

# The Tip of the Iceberg: How Pipette Tips Influence Results. Part 2: Perfect Geometry Makes a Difference

ULRIKE GAST, EPPENDORF AG, HAMBURG, GERMANY  
LAURA KOCH, EPPENDORF POLYMERE, HAMBURG, GERMANY

## Abstract

We performed a study including standard tips from different manufacturers in order to investigate the tip-related influences on the pipetting result. The study showed a dramatic impact of the tip on the pipetting accuracy. Small volume pipetting accuracy is significantly influenced by the quality of the tip orifice. Any impairment of geometry and shape can lead to retention of liquid.

## Introduction

Tip-related influencing factors on pipetting results can be compared to an iceberg. Some influencing factors are easy to observe, like the necessity for high tip attachment forces in order to achieve a secure and tight fit of non-recommended tips on the pipette cone. Other factors stay rather unknown, like the geometry of the tip orifice.

Part 1 of this series of Application Notes described pipette and tip as a system which may not perform within permissible error tolerances with all tips [1].

This is in concert with ISO 8655 [2] which regards pipette and tip as a system and the use of other tips than recommended as a system disruption. In part 1, the air-cushion size was detected to be the main influencing factor on bigger volumes like 1,000 µL. In contrast, with 10 µL tips other factors come into play. Here, the quality of the tip's orifice was found to be of high importance. For further investigation, the orifices of tips from different manufacturers were scrutinized and compared to calibration results [1,3].

## Material and methods

The orifice microscopy of 10 µL tips from Eppendorf and three other suppliers (manufacturers E, F, and H) was performed using a microscope (Leica®) with 25-fold magnification and a DFC 280 digital camera (Leica).

## Results and discussion

The zone where liquid leaves the tip during dispensing is very important for the accuracy of results. At this part of the tip, the drop cut-off occurs. Any imperfection of geometry or shape, e.g. by production errors, can lead to liquid retention. This especially plays a role with small volumes. A poor drop cut-off may not only impair the pipetting result but can make it impossible to dispense small volumes below 1 µL onto a solid surface.

Fig. 1 (see next page) shows examples of such production flaws. Tips from these manufacturers lead to the system performing outside permissible systematic error tolerances [1,3]. This underlines the direct influence of tip orifice geometry on pipetting of small volumes.

A good tip has a perfectly round geometry and uniform wall thickness. The orifice does not show flashes or loose material residues. Walls with varying thickness and not perfectly

round geometry lead to deflection of liquid drops to the outside of the thinner wall. Flashes cause liquid retention. Such production flaws are generated within tip production by tools with low quality or long maintenance cycles and/or a non-optimized injection molding process.

Since the production obviously has a high impact on the geometry and shape of the tip, influencing the pipetting result, it makes sense to take a closer look at production. The next part of this series will focus on the production of pipette tips and, in this context, the use of additives.

The thumbnail shows the article's title, authors, and abstract. It also features two scatter plots labeled '1 µL' and '10 µL'. The y-axis for both is 'Relative error (%)' and the x-axis is 'Systematic error (%)'. The 1 µL plot shows a wider spread of data points, while the 10 µL plot shows a more concentrated distribution. A red shaded area in both plots indicates the permissible error range.


## Missed part 1?

A rising number of published experiments cannot be reproduced by other groups. In general, only little attention is paid to the purchase of plastic laboratory consumables, like pipette tips and microtest tubes.

The fact that a pipette tip fits onto a pipette cone does not say anything about the pipetting accuracy of the pipetting system comprising the components "Pipette and Tip". We have studied this issue thoroughly in part 1 of our series "The Tip of the Iceberg: How Pipette Tips Influence Results".


If you missed part 1, then check out BioNews No. 44 in our archive at [www.eppendorf.com/bionews](http://www.eppendorf.com/bionews).

## The Tip of the Iceberg: How Pipette Tips Influence Results. Part 2: Perfect Geometry Makes a Difference



### Eppendorf epT.I.P.S.® 10 µL

The orifice has a good geometry and the function is not negatively influenced by production errors.




### Manufacturer E

**Problem 1:** Lying flashes caused by non-harmonized ejection molding process. Cavity has not been fully filled with liquid PP.

**Result:** Risk of deflection of water drop because of varying diameter of frontal area

**Problem 2:** Flashes at exterior wall


**Result:** Risk of keeping liquid residues



### Manufacturer F

**Problem 1:** Flashes at interior and exterior wall caused by long maintenance cycles of tool.

**Result:** Risk of liquid residues and risk of PP particles falling into sample. Additionally, particles inside the tip displace water leading to wrong liquid volume. Such an error pattern makes it impossible to dispense e.g. 0.8 µL sample onto a solid surface.



### Manufacturer H

**Problem 1:** Molding flashes caused by a non-tightening tool or problems with injection molding process where too much liquid PP has been injected.

**Result:** Liquid retention

**Problem 2:** Noticeable core shift. Not all walls have the same thickness. This error pattern can be caused by e.g. a poorly manufactured tool.

**Result:** Risk of deflection of liquid beside instead of into the target vessel.

**Fig. 1:** Microscopic pictures of 10 µL tip orifices of different manufacturers. Manufacturer E, F, and H failed the calibration by impaired systematic error. The examples have been chosen to explain production errors. PP=polypropylene.

### Literature

[1] The Tip of the Iceberg: How Pipette Tips Influence Results. Part 1: Tip Fit Is Not All Users Should Look for. BioNews 44, January 2016. [www.eppendorf.com/bionews](http://www.eppendorf.com/bionews)

[2] DIN EN ISO 8655:2002. Piston-operated volumetric apparatus. Beuth-Verlag, Berlin, Germany

[3] Application Note 354: The Tip of the Iceberg: How Pipette Tips Influence Results. [www.eppendorf.com/applications](http://www.eppendorf.com/applications)

### Readers' service

Eppendorf Liquid Handling Consumables • Ref. no. 288