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Literatur

- [1] Hussain M. 2016 'Molecular Imprinting' as Multidisciplinary Material Science: Today and Tomorrow. *UK J Pharmaceut Biosci*, 4: 121-132.
- [2] Campos J, Jiménez C, Trigo C, Ibarra P, Rana, Deepthi, Thiruganesh R, Ramalingam Murugan, Haidar Z S. 2015 Quartz Crystal Microbalance with Dissipation Monitoring: A Powerful Tool for BioNanoScience and Drug Discovery. *J Bionanosci*, 9: 249-260.
- [3] Hussain M, Northoff H, Gehring FK. 2015 QCM-D providing new horizon in the domain of sensitivity range and information for haemostasis of human plasma. *Biosens Bioelectron*, 66: 579-84. doi: 10.1016/j.bios.2014.12.003.
- [4] Dixon CM. 2008 Quartz Crystal Microbalance with Dissipation Monitoring: Enabling Real-Time Characterization of Biological Materials and Their Interactions. *J Biomol Tech*, 19: 151-158.
- [5] Sauerbrey G. 1959 Verwendung von Schwingquarzen zur Wägung dünner Schichten und zur Mikrowägung. *Z Phys*, 155: 206-222.
- [6] Kanazawa K, Gordon JG. 1985 The oscillation frequency of a quartz resonator in contact with a liquid. *Anal Chim Acta*, 175: 99-105.
- [7] Pomorska A, Shchukin D, Hammond R, Cooper MA, Grundmeier G, Johannsmann D. 2010 Positive frequency shifts observed upon adsorbing micron-sized solid objects to a quartz crystal microbalance from the liquid phase. *Anal Chem*, 82: 2237-2242.
- [8] Dybwad GL. 1985 A sensitive new method for the determination of adhesive bonding between a particle and a substrate. *J Appl Phys*, 58: 2789-2790.
- [9] Gehring, FK. 2006 Schwingquarzensensorik in Flüssigkeiten -Entwicklung eines Blutanalysegerätes, Cuvillier Verlag, Göttingen, ISBN 9783865378729.
- [10] Costerton JW, Stewart PS, Greenberg EP. 1999 Bacterial biofilms: A common cause of persistent infections. *Science*, 284: 1318-1322.
- [11] Singh PK, Schaefer AL, Parsek MR, Moninger TO, Welsh MJ, Greenberg EP. 2000 Quorum-sensing signals indicate that cystic fibrosis lungs are infected with bacterial biofilms. *Nature*, 407: 762-764.
- [12] Boucher HW, Talbot GH, Bradley JS, Edwards JE, Gilbert D, Rice LB, Scheld M, Spellberg B, Bartlett J. 2009 Bad Bugs, No Drugs: No ESKAPE! An Update from the Infectious Diseases Society of America. *Clin Infect Diseases*, 48: 1-12.
- [13] Olsson AL, Mitzel MR, Tufenkji N. 2015 QCM-D for non-destructive real-time assessment of *Pseudomonas aeruginosa* biofilm attachment to the substratum during biofilm growth. *Colloids Surf B Biointerfaces*, 136: 928-934. doi:10.1016/j.colsurfb.2015.10.032.
- [14] Olofsson AC, Hermansson M, Elwing H. 2005 Use of a quartz crystal microbalance to investigate the antiadhesive potential of Nacetyl-L-cysteine. *Appl Environ Microbiol*, 71: 2705-2712.
- [15] Olsson ALJ, van der Mei, HC, Busscher HJ, Sharma PK. 2008 Influence of cell surface appendages on the bacterium-substratum interface measured real-time using QCM-D. *Langmuir* 25: 1627-1632.
- [16] Sismaet HJ, Abadian PN, Goluch ED. 2014 Monitoring Bacterial Biofilm Growth and Removal. Department of Chemical Engineering, Northeastern University, Boston, MA, USA. *Appl. Note*.
- [17] Castelain M, Koutris, E, Andersson M, Wiklund, K, Björnham O, Schedin, S, Axner O. 2009 Characterization of the Biomechanical Properties of T4 Pili Expressed by *Streptococcus pneumoniae* – A Comparison between Helix-like and Open Coil-like Pili. *Chem Phys Chem*, 10: 1533-1540, DOI: 10.1002/cphc.200900195.
- [18] Belas R, 2014 Biofilms, flagella, and mechano-sensing of surfaces by bacteria. *Trends Microbiol* 22: 517-527.
- [19] Zhang S, 2003 Fabrication of novel biomaterials through molecular self-assembly. *Nat. Biotechnol*, 21, 1171-1178.
- [20] Jessel N, Oulad-Abdelghani M, Meyer F, Lavalle P, Haykel Y, Schaaf P, Voegel JC. 2006, *Proc Natl Acad Sci US A*, 103, 8618-8621.
- [21] Waslathanthri DP, Kuhn L, Rusling JF. 2013 Multilayer Bio-Films for bone tissue regeneration. *University of Connecticut Health Center*, Storrs CT. *Appl. Note*.
- [22] Eap S, Richert L, Lemoine S, Kalaskar D, Demoustier-Champagne S, Atmani H, Mély Y, Fioretti F, Schlatter G, Kuhn L, Ladam G, Benkirane-Jessel N. 2014 Osteogenetic properties of electrospun nanofibrous PCL scaffolds equipped with chitosan-based nanoreservoirs of growth factors. *Macromol Biosci*. 14: 45-55. doi: 10.1002/mabi.201300283.